



POLYNOR AS has a long and proud history. Since the founder Odd Terje Østgaard started POLYNOR AS in 1994, the company has gained international recognition for its patented product POLYSAFE® Safety Box and Incineration Container.

POLYNOR's range of safety boxes has since 1996 been tested and approved by an independent WHO-accredited laboratory to comply with the requirements of WHO Standards E10/IC.1 and E10/IC.2 and WHO/PQS/E10/SB01.1 on needle penetration, water resistance and combustibility etc..

Today the company consists of highly motivated and service minded employees with long experience. They are known for their efficiency and quick service.



POLYNOR AS · Brennerigt 3 · PO Box 1273 · N-2806 Gjøvik / Norway
 Tel: +47 61 13 89 30 · Fax: +47 61 13 89 40
 Email: post@polynor.no · www.polynor.no
 Org.nr: NO 968 729 012 MVA

A GREEN COMPANY


POLYNOR AS is an environmentally friendly company. All POLYSAFE® Safety Boxes are produced from recycled solid board. They are flat packed which saves on transportation costs and fuel consumption.

FAST DELIVERIES

The boxes are produced in Norway in a high-tech plant. ISO Certified. We provide fast and secure production and delivery all over the world.

SAFE INJECTIONS

POLYSAFE® Safety Boxes prevent health workers and communities from being exposed to unnecessary risks from contaminated needles.



POLYNOR
 The maker of POLYSAFE®
 Safety Boxes since 1994

POLYSAFE® Safety Boxes for safe and easy collection and destruction of used syringes and needles



POLYSAFE® Safety Boxes are approved by WHO

The maker of
POLYSAFE®
Safety Boxes
since 1994



**POLYSAFE® Safety Boxes
are well-known products
in the emergency relief
market all over the world**

Among our clients are:

- ✓ Major international health agencies
- ✓ Major regional health agencies
- ✓ Non-governmental organizations (NGO'S)
- ✓ National health ministries
- ✓ US Military (NSN 6530-01-460-4782)



**Safety boxes have many
areas of use**

- ✓ Vaccination campaigns
- ✓ Hospitals
- ✓ Mobile health clinics
- ✓ Medic in the armed forces

**In the fight to save lives,
use POLYSAFE® Safety
Boxes to collect and
destroy used syringes and
needles**



**A POLYSAFE® Safety Box
which is in itself an
incinerator and the only
one to be burned on-site**

On site destruction

- ✓ Use under primitive conditions in which the objective is to get the used syringes as fast as possible out of circulation.
- ✓ Can be used anywhere
- ✓ Unnecessary with transportation to a suitable incineration site
- ✓ Does not require lighter fluid
- ✓ National Stock Number (NSN) 6530-01-460-4782



POLYSAFE® 5L Safety Box & Sharps Container

WHO CODE PQS E10/01,
printing in English, French, Spanish and Russian, flat packed in robust water resistant shipping cartons of 25 each.



POLYSAFE® 10L Safety Box & Sharps Container

WHO CODE PQS E10/02,
printing in English, French, Russian and Spanish, flat packed in robust water resistant shipping cartons of 15 each.



POLYSAFE® 15L Safety Box & Sharps Container

WHO CODE PQS E10/03,
printing in English, French, Russian and Arabic, flat packed in robust water resistant shipping cartons of 10 each



POLYSAFE® 20L Safety Box & Sharps Container

WHO PIS CODE E12/09,
approved according to WHO/UNICEF Standard E10/IC.2
printing in English, French, Spanish, Russian and Arabic, flat packed in robust water resistant shipping cartons of 10 each.



POLYSAFE® 2,5L Safety Box & Sharps Container

WHO PIS CODE E12/17,
approved according to WHO/UNICEF Standard E10/IC.2
printing in English and Russian, flat packed in robust water resistant shipping cartons of 30 each.



POLYSAFE® 5L Safety Box & Incineration Container

WHO PIS CODE E12/02,
approved according to WHO/UNICEF Standard E10/IC.1 and E10/IC.2,
NSN 6530-01-460-4782
printing in English, French, Russian and Spanish, flat packed in robust water resistant shipping cartons of 25 each.

POLYNOR is
the only
manufacturer
of multiple
sizes of
Safety Boxes