




POLYSAFE® 5L Safety Box and Incineration Container



POLYNOR

The maker of
POLYSAFE® Safety Boxes since 1994

Safety box for the disposal of used syringes and needles WHO CODE PIS E12/02

SPECIFICATIONS

Volume	5 litres
Capacity	155 autodis. syringes of 0.5 ml with needles and caps
External dimensions:	
Before assembling	590x272x3mm
After assembling	285x160x130mm
Materials	recycled board, waterproof
Thickness of walls	1.4 mm
Weight fully assembled	310.0 grs
Diam. syringe insert whole	35 mm
Colour	White box with black and red print

COMMENTS

*Test Reports: Force Institute 65788/m1460-146 (1996)
and 91530/m1460-190.c (1998) and 128258/ml1460- 284.b (2000)
Meets WHO/UNICEF Standard E10/IC.1.*

Directions for use on the box in English, French, Spanish and Russian. Other languages on request.
Boxes of 25 pcs. flat packed in robust water resistant outer carton.

Shipping volume 0.02 m³

Shipping weight (25 boxes) 9.0 kg

POLYNOR AS · Brennerigt 3 · PO Box 1273 · N-2806 Gjøvik / Norway

Tel: +47 61 13 89 30 · Fax: +47 61 13 89 40 · Email: post@polynor.no · www.polynor.no

Foretaksregisteret NO 968 729 012 MVA