




POLYSAFE® 20L Safety Box and Sharps Container



POLYNOR

The maker of
POLYSAFE® Safety Boxes since 1994

Safety box for the disposal of used syringes and needles WHO CODE PIS E12/09

SPECIFICATIONS

Volume	20 litres
Capacity	654 autodis. syringes of 0.5 ml with uncapped needles
External dimensions	
Before assembling	946 x 463 x 2.8mm
After assembling	405 x 285 x 193mm
Materials	recycled board
Thickness of walls	1.4 mm
Weight fully assembled	840.0 g
Diam. syringe insert hole	35 mm
Color	White box with black and red print

COMMENTS

Test Reports: Force Institute 128258/m1460-284.d(2000)

Meets WHO/UNICEF Standard E10/IC.1.

Directions for use on the box in English, French, Spanish and Russian, or other languages on request.

Boxes of 10 pcs. flat packed in robust water resistant outer carton.

Shipping volume 0.02 m³

Shipping weight 8,1kg

POLYNOR AS · Brennerigt 3· PO Box 1273 · N-2806 Gjøvik / Norway

Tel: +47 61 13 89 30 · Fax: +47 61 13 89 40 · Email: post@polynor.no · www.polynor.no

Foretaksregisteret NO 968 729 012 MVA